



CK SURGITECH

FACE MASK LEVEL 3

4 ply - DISPOSABLE

TGA
APPROVED


50 masks / box
Order Code : **CK-M001**

INNER LAYER: SOFT AND WEARABLE FOR USER COMFORT

2ND LAYER: MINIMISES PARTICLE DISPERSAL

3RD LAYER: DEFLECTS FLUID

OUTER LAYER: NON-WOVEN BREATHABLE



Characteristic

Level 3

Bacterial Filtration Efficiency	>98%
Sub-micron particulates filtration efficient at 0.1 micron	>98%
Differential pressure, mm H2O/cm2 breathability	<5.0
Resistance to penetration by synthetic blood (fluid resistance) min pressure in mmHg for pass result	160mm HG

CK Surgitech Pty Ltd

6/14 Technology Drive, Arundel
Queensland Australia 4214.
P: 07 5574 5573 E: sales@cksurgitech.com