



CK SURGITECH

OC2 - PACK

Olympus® cleaning duo pack

Order Code :
GAR -03875NS

25 pcs/ box

andorate®



Suction Channel
Cleaning Adapter



Air Water Cleaning Adapter

MEET YOUR SINGLE USE INFECTION CONTROL NEEDS IN ONE CONVENIENT PACK

Cleaning is a critical step in endoscope reprocessing that starts in the procedure room. The OC2-PACK by Andorate™ contains both the Air Water Channel Cleaning Adapter and the Suction Channel Cleaning Adapter for Olympus® GI Endoscopes. The Air Water Channel Cleaning Adapter flushes the endoscopes air channel with water, preparing it for reprocessing. The suction channel cleaning adapter is used to aspirate reprocessing fluids from the distal end of the Endoscope through the instrument channel including the biopsy port. The Andorate™ OC2-PACK is a disposable alternative to the reusable options.

Features

- Eliminates the need to track, trace and reprocess post-procedure
- Complies with infection control guidelines
- One convenient pack

Olympus® is the registered trademark of Olympus Corporation

CK Surgitech Pty Ltd

6/14 Technology Drive, Arundel
Queensland Australia 4214
P: 07 5574 5573 E: sales@cksurgitech.com