



CK SURGITECH

FACE MASK - 4 ply DISPOSABLE

50 masks / box

Order code **CK-M001/P**




INNER LAYER: SOFT AND WEARABLE FOR USER COMFORT

2ND LAYER: MINIMISES PARTICLE DISPERSAL

3RD LAYER: DEFLECTS FLUID

OUTER LAYER: NON-WOVEN BREATHABLE



Characteristic

Specification

| | |
|--|----------|
| Bacterial Filtration Efficiency | 99% |
| Differential Pressure (Delta P) mm H ₂ O/cm | <5.0 |
| Fluid Resistance Penetration | 160mm HG |
| Particulate Filtration Efficiency | 99% |

CK Surgitech Pty Ltd

6/14 Technology Drive, Arundel
Queensland Australia 4214
P: 07 5574 5573 E: sales@cksurgitech.com