



CK SURGITECH

# COLD SNAP Polypectomy Snare

Order Code :  
**AG-5336-241023**

**10** pcs/ box



**FOR THE ULTIMATE IN CONTROL AND SUPERIOR CUTTING**

## Features

## Potential Advantages

Cold Hexagonal 10mm Snare	The hexagonal shape allows for the manipulation and control of tissue and polyps. 10mm Loop is effective for the resection of small and diminutive polyps. Cold polypectomy enables a superior cut with cleaner margins.
Firm 0.27mm Wire	Designed for users who prefer a firmer wire loop for cold polypectomy.
Ergonomic Handle	For ease of use and user comfort
Full 1 to 1 Rotation	Ultimate control, facilitating decreased procedure time which can support decreased anaesthetic time

**CK Surgitech Pty Ltd**

6/14 Technology Drive, Arundel  
Queensland Australia 4214  
P: 07 5574 5573 E: sales@cksurgitech.com