



CK SURGITECH

DUOMASTER Polypectomy Snare


Order Code :

AG-5076-241023 - 10mm Snare

AG-5076-241523 - 15mm Snare



10 pcs/ box



FOR THE ULTIMATE CONTROL HOT OR COLD

Features

Potential Advantages

Hexagonal 10mm & 15mm Snare	The hexagonal shape allows for the manipulation and control of tissue and polyps.
Hybrid	Designed to perform optimally for both hot and cold snare polypectomy.
Thin/ stiff 0.3mm wire 1 x7 formation	Facilitates a superior clean cut for cold polypectomy. Maintains loop shape even after multiple passes.
Ergonomic Handle	For ease of use and user comfort
Full 1 to 1 Rotation	Ultimate control, facilitating decreased procedure time which can support decreased anaesthetic time

CK Surgitech Pty Ltd

6/14 Technology Drive, Arundel
Queensland Australia 4214

P: 07 5574 5573 E: sales@cksurgitech.com