




CK SURGITECH

VARIO HOT Polypectomy Snare

Order Code :
PFS02-33023230

10 pcs/ box

THE ULTIMATE 3 in 1 HOT SNARE



Features

Benefits

Combination Hot Snare 10/20/30mm	<ul style="list-style-type: none"> • Eliminates the need to exchange snares during a procedure. • Easily manipulated between small and larger sized polyps, reducing procedure time and cost of multiple snares.
Stiff 0.4mm wire 1x7 formation	<ul style="list-style-type: none"> • The loop design maintains its shape and intensity after multiple resections. • The stiffer wire minimises distal deflection.

CK Surgitech Pty Ltd

6/14 Technology Drive, Arundel
Queensland Australia 4214
P: 07 5574 5573 E: sales@cksurgitech.com