



CK SURGITECH

DIAMOND CUT COLD SNARE

10 pcs/ box

**SUPERIOR CUT
CLEANER MARGINS**



CK Surgitech Pty Ltd

6/14 Technology Drive, Arundel
Queensland Australia 4214
P: 07 5574 5573 E: sales@cksurgitech.com



DIAMOND CUT COLD SNARE

| Product Code | Description | Loop size | Length (cm) | Sheath Diameter (mm) |
|--------------|-------------------|-----------|-------------|----------------------|
| CS3-11023230 | Snare Diamond Cut | 10mm | 230 | 2.3 |
| CS3-11523230 | Snare Diamond Cut | 15mm | 230 | 2.3 |

Features

Benefits

| | |
|-------------------------|---|
| Dedicated Cold Snare | <ul style="list-style-type: none">• Specifically designed for resecting small and diminutive polyps using a cold snare technique• Cold polypectomy enables a superior cut with cleaner margins |
| Diamond Shape 10mm/15mm | <ul style="list-style-type: none">• The diamond shape loop and flexible wire captures smaller polyps with ease |
| Thin wire | <ul style="list-style-type: none">• The thin 0.23mm cutting wire provides a guillotine cut that reduces the force required to complete the resection and that may limit tearing or snagging |
| Ergonomic Handle | <ul style="list-style-type: none">• Accommodates various hand sizes• The three-ring handle features a shortened throw and comfortable grip |
| Stiff Catheter | <ul style="list-style-type: none">• The stiff catheter aids in insertion through an extended scope position while providing support for a smooth effortless cut |