



CK SURGITECH

# TECHBITE BIOPSY FORCEPS

Standard Capacity

10<sub>pcs</sub>/ box

SHARPER CUP  
IMPROVED BITE  
OPTIMAL CONTROL



CK Surgitech Pty Ltd

6/14 Technology Drive, Arundel  
Queensland Australia 4214  
P: 07 5574 5573 E: sales@cksurgitech.com



CK SURGITECH

# TECHBITE BIOPSY FORCEPS

## Standard Capacity

Product Code	Description	Length (cm)	Sheath Diameter (mm)
EBF33-11023230	Biopsy forceps TECH BITE oval cup plain	230	2.3
EBF33-11123230	Biopsy forceps TECH BITE oval cup plain with spike	230	2.3
EBF43-11023230	Biopsy forceps TECH BITE oval cup serrated	230	2.3
EBF43-11123230	Biopsy forceps TECH BITE oval cup serrated with spike	230	2.3

### Features

### Benefits

Large spoon cup design	<ul style="list-style-type: none"><li>• Quality sample</li><li>• Superb cutting performance</li></ul>
Coated Sheath	<ul style="list-style-type: none"><li>• Reduces sliding friction</li><li>• Protects the endoscope working channel</li></ul>
Thumb ring / short throw handle	<ul style="list-style-type: none"><li>• Provides optimal control</li><li>• 1 to 1 response tactility</li></ul>
Range of cups- plain and serrated with and without spike	<ul style="list-style-type: none"><li>• Meets all clinical requirements</li></ul>



Spike



Serrated



Plain

CK Surgitech Pty Ltd

6/14 Technology Drive, Arundel  
Queensland Australia 4214  
P: 07 5574 5573 E: sales@cksurgitech.com