


LOCK ONTO TISSUE **LOCKADO™ HEMOSTASIS CLIP**


Hemostasis clips need to provide accuracy, reliability and precise functionality. Lockado allows unrestricted opening and closing so it can be repositioned until the desired position is achieved.

With two sets of opposing atraumatic teeth in the arms of the clip, Lockado delivers a reliable and secure grip for consistent performance. Driven by Micro-Tech Endoscopy's proven clip delivery system, the Lockado clip provides exceptional reliability and tissue management.

- Jaw style works for multiple clinical applications
- Can be used with side-viewing endoscope
- Shorter stem make clips less obtrusive
- Variety of sizes: 11 mm, 16 mm and **THE NEW 20 mm and 8mm** arm opening widths



KEY FEATURES



SPECIFICATIONS

LOCKADO HEMOSTASIS CLIP, REPOSITIONABLE, ROTATABLE

| Order Number | Opening Width (mm) | Sheath Diameter (mm) | Working Length (cm) | Minimum Channel Size (mm) | Packing Units |
|---------------|--------------------|----------------------|---------------------|---------------------------|---------------|
| LOCK-C-26-230 | 8 | 2.6 Pink | 235 | 2.8 | 10/Box |
| LOCK-D-26-230 | 11 | 2.6 Blue | 235 | 2.8 | 10/Box |
| LOCK-F-26-230 | 16 | 2.6 Yellow | 235 | 2.8 | 10/Box |
| LOCK-G-26-230 | 20 | 2.6 Grey | 235 | 2.8 | 10/Box |